


# SLEEPONE PROSV

Cheyne-Stokes Respiration, Central Apnea, Obstructive Apnea

## All-in-One Solution


**Servo Ventilation + Auto Cpap + Respiratory Backup = SleepOne ProSV™**

- Auto Pressure Support
- Auto EPAP
- Auto Backup Rate
- Can also work as CPAP Mode Servo Ventilation (Min Pressure Support = 0)



Cheyne-Stokes Respiration

A Complex Device for increased Patient Complexity !



ProSV gives instant reaction to changes in patient effort.

### Technical Specifications

Pressure Range	4-35 CmH <sub>2</sub> O
Working Mode	ProSV Active Servo Ventilation, CPAP
Ventilation Type	Leakage Ventilation
Voltage	100 - 240 V AC
Humidifier	SleepOne Heated Humidifier
Data management	SD Card stores up to 365 nights of detailed therapy data.



### Event

Obstructive Apnea  
Hypopnea  
Vibratory Snore  
Patient Flow Decrease  
Complex Apnea  
RERA  
High Leak  
Periodic Breathing  
Low Respiratory Rate

### Action

EPAP Change  
EPAP Change  
EPAP Change  
Pressure Support Change  
Airway blockage test  
EPAP Change  
Marked for Reporting  
Marked for Reporting  
Backup Breath



# Active Servo Ventilation Device



## **Active Patient Control on Breath-by-breath basis**

*SleepOne ProSV™ is an active and also proactive device which responds to patients breath by breath changes. Patient who suffers from Cheyne-Stokes Ventilation, Central Apnea, Obstructive Apnea or Complex Apneas can benefit from the active/proactive algorithm of ProSV™. Device monitors several parameters from the patients and adjust pressure and breath rate according to the change. Patients who don't get enough from standard Bilevel therapy, or Auto-Bilevel algorithms can find the ProSV™ is a problem solver for complex patients.*

## **Patient Focused Flexible Algorithm**

*ProSV™ targets to keep the patient flow to the %90 of the 3 minute moving average. During a normal breathing cycle, in which patient's respiratory effort is minimum, device keeps the pressure level to a minimum also. If there is no respiratory effort, patient is forced with the backup rate. According to the measured patient flow, ProSV™ categorizes the events. During complex apneas, device does not respond. Obstructive apneas are responded with a pressure increase after multiple events. Hypopneas are also responder events with pressure increase after multiple events. Snores are also responder events. Pressure increase is monitored and the amount is interactively chosen by the algorithm. ProSV™ decreases the pressures during no-event times. Patient is always treated with necessary pressure. High Leak Flows, Periodic breathing (Cheyne-Stokes Ventilation) and all other events are recorded for reporting.*