





OXYGEN CONCENTRATORS

Kare Medical proudly presents the new Oxybreath Mini series, Mini 3 and Mini 5.

Both with the non medical look, low volume design that saves transportation costs of providers and makes the therapy easier for patients. Mini 3 is specially designed for low flow patients in order to save them from the cost of 5 liter machines. Mini 5 is for patients that require regular therapy in a smaller and power efficient platform.



Part Number
Operating Voltage
Mains Frequency
Power
Oxygen Flow
Oxygen Saturation

O₂ Monitor

Class

OXY6300 220-240V / AC 50-60 Hz. 165 W 0,5-3 LPM 1-3 L 93 ± 3%

- High Oxygen %O₂ > 82% - Low Oxygen %O₂ < 82%

- Technical Service Required < 70%

Operating Pressure Operating Temp. Gross Weight Net Weight Type

0,07-0,09 Mpa 5-40° C 14,4 Kgs 12 Kgs. B

II





3-5











Key Points

- Ultra low volume decreases shipping costs (288 pcs in a 20' container)
- Very small visual appearance and non medical look
- Low power consumption&less heat production
- 3 liter version is available for low flow patients, saves from unnecessary power cost and extra weight of 5 liter machines
- No maintenance required for 2 years
- Silent operation

Mini 5 Technical Specifications

Part Number OXY6500
Operating Voltage 220-240V / AC
Mains Frequency 50-60 Hz.
Power 330 W
Oxygen Flow 0,5-5 LPM
Oxygen Saturation 1-5 L 93 ± 3%

 O_2 Monitor - High Oxygen $\%O_2 > 82\%$ - Low Oxygen $\%O_2 < 82\%$

- Technical Service Required < 70%

Operating Pressure 0,07-0,09 Mpa
Operating Temp. 5-40° C
Gross Weight 17,5 Kgs
Net Weight 15,1 Kgs.
Type B
Class II

