

AEROCARE II

PISTON TYPE COMPRESSOR NEBULIZER

Hikoneb® AeroCare II Piston Type Nebulizer is designed to be used on living environment of the patient. AeroCare II unit is perfect for daily usage of a single patient.

Nebulizer is a device used to administer medication to people in forms of a liquid mist to the airways. It is thus used in treating respiratory diseases. Piston type nebulizers pump air through a liquid medicine to turn into a vapor, which is then inhaled by the patient.

These medicines are frequently steroids and the reason they are inhaled instead of digested is to limit their effect to the lungs / respiratory system only. Otherwise, that amount of steroid saturated throughout the body would be toxic.

Technical Specifications

Operating Voltage	230V / AC
Mains Frequency	50 Hz.
Power	70 W
Maximum Pressure	3 bar
Nebulization Rate	0-5 ml/min
Weight	1.75 Kg
Type	B
Class	II



Main Areas of Usage

- Drug Nebulization
- Specific - Non specific lung diseases
- Bronchial Asthma
- Pneumonias
- Acute and Non-acute bronchial

