

Aerocare Pediatric

A101- A102 PISTON TYPE NEBULIZER

Hikoneb® AeroCare 101 and 102 Piston Type Nebulizer is designed to be used for pediatric patients. Both of the units have attractive design that makes the therapy more enjoyable for kids.



A101



A102

Main Areas of Usage

- Drug Nebulization
- Specific - Non specific lung diseases
- Bronchial Asthma
- Pneumonias
- Acute and Non-acute bronchial

Technical Specifications

Operating Voltage	230V / AC
Mains Frequency	50Hz
Power	70 W
Maximum Pressure	3 bar
Nebulization Rate	0-5 ml/min
Weight	(101) 1.42 KG (102) 1.4 KG
Type	B
Class	II

