


# ECOASPIR

## PORTABLE SUCTION UNITS

ECOASPIR series compact suction units specially designed to be lightweight and comfortable for home use. ECOASPIR has 1000 ml. disposable collection bottle with automatic float shut-off. It meets the ISO 10079 standards.

### Main Features

- High Vacuum Capability, 550 mmHg
- Elegant, Small and Light weight, Portable Design
- Integrated Bacterial Filter
- Automatic Float Shut-off
- 1000 ml. Disposable Collection Bottle
- Available both battery operated and mains operated models



### Mains Operated Model Technical Specifications

Weight	2,200 gr. *
Dimensions	240 x 270 x 145 mm ( W x D x H )
Collection Bottle	W x D x H )
Volume	1000 ml.
Float Safety	Integrated Float Valve
Max. Vacuum	550 mmHg **
Min. Vacuum	100 mmHg **
Max. Flow	9 Liter/min **
Voltage	220 V/AC 50Hz
Power	80 Watt

### Battery Operated Model Technical Specifications

Weight	3000 gr. *
Dimensions	240 x 270 x 145 mm ( W x D x H )
Collection Bottle	W x D x H )
Volume	1000 ml.
Float Safety	Integrated Float Valve
Max. Vacuum	590 mmHg **
Min. Vacuum	100 mmHg **
Max. Flow	10 Liter/min **
Voltage	230 V/AC 50-60 Hz
Power	45 Watt