

SLEEPONE BILEVEL ST AUTO



All in One Solution for Developing Markets...



Modes Included

- 1- Cpap
- 2- AutoCpap
- 3- Bilevel S
- 4- Bilevel T
- 5- Bilevel ST
- 6- Bilevel Auto Fixed PS
- 7- Bilevel Auto Variable PS
- 8- Bilevel Auto ST Fixed PS
- 9- Bilevel Auto ST Variable PS

SleepOne Bilevel ST Auto is the newcomer from the SleepOne Product Line. With its powerful blower device can go up to 35 cmH₂O in allowed modes.

It is a solution developed for developing countries. Single product fits wide range of patients. This aspect makes the product charming to the professional users. When the patient is ready to go home, or a less complicated solution is needed after the intensive care, SleepOne Bilevel ST Auto provides the freedom to try different modes on the patient to choose which one is the best.



*Together with the SleepOne Desktop software, all the settings can be changed, modes can be changed or the patient can be monitored realtime with several graphics**

Graph 1 : Pressure

Graph 2 : Flow Graph including Patient Flow Curve, Leak Flow Curve, Event Notifiers (Apnea, Hypopnea, etc.)

Graph 3: Lung Fill Curve

All the curves give you realtime data at your fingertips. No additional hardware needed, other than a USB port and a Windows PC.

Data management : SD Card stores up to 365 nights of detailed therapy data.



• Realtime monitoring and Graphics are done via SleepOne Desktop Software

