

SLEEPONE HEATED HUMIDIFIER



SleepOne Humidifier helps for preventing the symptoms of dryness that most Cpap users experience. User friendly designed truly integrated unit allows ultimate patient compliance and increases the chance of therapy success.



Key Features

- Easy to clean
- Easy to adjust via Cpap Menu, even in the dark
- Compact design, easily fits to SleepOne Cpap bag
- High volume chamber

Technical Specifications

Dimensions	262x159x130 mm (LXWXH)
Weight	51 gr
Chamber Volume	600 ml
Chamber	ABS Lid, Polycarbonate Body, Stainless Steel Plate (SS 316)
Docking Station	Thermoplastic
Heater Plate	Stainless Steel(SS 316)
Voltage	100-230V AC,50/60Hz
Power	85W
Heater Plate Temperature Range	30-65° C
Heater Plate Safety Cutout	105° C

